

# CAUTION: Choosing Your Help Button

It is important that you and your caregiver choose the Personal Help Button most suitable for you. Please consider the following:

## **Pendant styles**

The Classic and Slimline Pendants are worn around the neck; this helps make the button accessible by either of your hands.

**NOTE: the Pendant's neck cord is not designed to break away. Therefore, it can pose a choking risk,** particularly to wearers in wheelchairs, using walkers, using beds with guard rails, or otherwise near protruding objects upon which the cord can become tangled. Wearers for whom this is a concern may wish to consider the Wrist style button.

## **Wrist styles**

The Slimline Wristband and Tempo Watch are worn around your wrist. In circumstances where the user's arm that is **not** wearing the pendant becomes immobilized (examples: due to stroke or arm being trapped under a heavy object), this can present a risk that the user will be unable to press the button.

**Please consult your caregiver or Lifeline Program on which Help Button style is most appropriate for you.**

**PHILIPS**  
**Lifeline**